

Polyclonal Antibody to Equine IgG (T) -FITC

Alternate names:	Horse IgG, Horse Immunoglobulin G
Catalog No.:	SP2007F
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host:	Goat
Immunogen:	Purified equine IgG (T).
Format:	State: Liquid purified IgG fraction. Purification: Affinity Chromatography. Buffer System: PBS, pH 7.2 containing 0.09% Sodium Azide as preservative and 0.2% BSA as stabilizer. Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry: 1/20-1/100. Immunohistochemistry on Frozen Sections: 1/20-1/100. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG (T) and shows no cross - reactivity with other immunoglobulin classes. Species: Horse. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com