

Polyclonal Antibody to Feline IgG -FITC

Alternate names:	Cat IgG, Immunoglobulin G
Catalog No.:	SP2013F
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host:	Goat
Immunogen:	Purified Feline IgG Fc fragment.
Format:	State: Liquid purified IgG fraction Purification: Affinity Chromatography Buffer System: PBS containing 0.09% Sodium Azide as preservative and 0.2% BSA as stabilizer Label: FITC – Fluorescein Isothiocyanate Isomer 1
Applications:	Flow Cytometry. Immunohistochemistry on Frozen Sections: 1/200-1/2,000. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises Feline IgG (Fc) and shows no cross-reactivity with other Feline immunoglobulin classes in immunoelectrophoresis. Species: Cat. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.
General References:	1. De Craeye, S. et al. (2008) Prevalence of Toxoplasma gondii infection in Belgian house cats. Vet Parasitol. 157:128-32. 2. Saari, S. et al. (2008) A nonhealing ulcerative skin lesion associated with Trichinella nativa infection in a cat. J Vet Diagn Invest. 20:839-43.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com