

Polyclonal Antibody to Chicken IgG (Fc) -Purified

Alternate names:	Chicken IgY
Catalog No.:	SP2021P
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host:	Goat
Immunogen:	Purified chicken IgG - Fc fragment
Format:	State: Liquid purified IgG fraction. Purification: Affinity Chromatography. Buffer System: PBS, pH 7.2 containing 0.09% Sodium Azide as preservative.
Applications:	ELISA (1/100-1/10000). Immunodiffusion. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises IgG(Fc) and shows no cross-reactivity with other chicken immunoglobulin classes in Immunoelectrophoresis. Species: Chicken. Other species not tested.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com