

Polyclonal Antibody to Porcine IgG (H&L) - HRP

Alternate names:	Pig IgG, Pig Immunoglobulin G
Catalog No.:	SP2086HRP
Quantity:	1 mg
Concentration:	1.0 mg/ml
Host:	Goat
Format:	State: Liquid purified IgG Purification: Affinity chromatography Buffer System: PBS containing 0.09% Proclin (TM)300 as preservative and 0.2% Bovine Serum Albumin as stabilizer Label: HRP – Horseradish Peroxidase
Applications:	ELISA: 1:10,000 - 1:100,000. Other applications not tested. Optimal dilutions are dependent on conditions and should be determined by the user.
Specificity:	This antibody recognises porcine IgG heavy and light chains and has been shown to react with light chains common to other pig immunoglobulins by immunoelectrophoresis and ELISA. The antibody may cross react with IgG from other species.
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com