



AVISCERA BIOSCIENCE

Rabbit Anti Secreted Frizzled-Related Protein 5 (SFRP5) (Human) IgG

Product Information

Code	A00179-01-100
Name	Rabbit Anti SFRP5 (H) IgG
Clone No.	Polyclonal
Lot No.	
Size	100 µg
Species	Human
Host	Rabbit
Immunogen	Human SFRP5
Ab Type	IgG
Purification	Protein A
Formulation	Lyophilized Form without Preservatives
Carry	Free
Storage	-20 ° C
Specificity	Human
Reconstitution	PBS, 100 µl
Application	IHC ELISA

PREPARATION

This antibody was produced from a rabbit immunized with purified, *E. coli*-derived, recombinant human SFRP5. That IgG was purified by Protein A affinity

FORMULATION

100 µg of Rabbit Anti Human SFRP5 purified IgG in 100 µl of PBS without preservatives was lyophilized.

RECONSTITUTION

Add 100 µl of PBS to the vial to prepare antibody stock solution at 100 µg/100 µl. Store reconstituted antibody at 2 to 8 ° C for up a few weeks. This antibody can also be aliquotted (by 10 µL per vial) and stored frozen at -20° C to -70° C **in a manual defrost freezer** for up six months without detectable loss of activity.

STORAGE

Lyophilized antibody can be stored at 2 ~8 ° C for a few weeks or at -20 ~ -70 ° C for 12 months. **Avoid repeated freeze-thaw cycles.**

SPECIFICITY

This antibody has been selected for its ability to recognize human SFRP5 in indirect ELISAS and immunohistochemistry.

APPLICATIONS

Indirect ELISA - This antibody can be used at 1:4000 (0.25 µg/ml) to detect recombinant Human SFRP5 on indirect ELISA.

Immunohistochemistry- That antibody can be used at 2-4 µg/ml to detect SFRP5 in human pancreas cancer tissues slides.

Optimal dilutions should be determined by each laboratory for each application.