




## Product Data Sheet

<b>Product Name:</b>	Anti-MMP-7 (CT)
<b>Catalog Number:</b>	29577
<b>Lot Number:</b>	See label on vial
<b>Product Description:</b>	This polyclonal antibody is supplied as an epitope affinity purified rabbit IgG, 50 µg in 250 µl of 1X PBS (pH 7.4) containing 0.1% BSA and 0.05% sodium azide.
<b>Immunogen:</b>	Rabbit anti-MMP-7 polyclonal antibody was raised against a synthetic peptide corresponding to the C-terminus of human MMP-7.
<b>Species Reactivity:</b>	Species reactivity includes human, mouse and rat, while others remain unknown.
<b>Application Notes:</b>	The following concentration ranges are recommended starting points for this product. Optimal working concentrations should be determined by the investigator for specific applications.
	WB: 0.5-2.0 µg/ml
	IHC: 10.0-20.0 µg/ml



This product is for *in vitro* research use only.