



Product Data Sheet

Product Name: β -Amyloid (3-40)

Catalog Number: AS-61029 (0.1 mg) Lot Number: See label on vial
AS-61029-1 (1.0 mg)

Sequence: H-Glu-Phe-Arg-His-Asp-Ser-Gly-Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly-Val-Val-OH (3-letter code)
EFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVV (1-letter code)

Molecular Weight: 4143.7

Peptide Purity: >95%

Appearance: Lyophilized white powder

Peptide Reconstitution: β -Amyloid (3-40) peptide is freely soluble in 1% NH₄OH.

Storage: β -Amyloid (3-40) peptide is shipped at ambient temperature. Upon receipt, store lyophilized peptide at -20°C or lower. Reconstituted peptide can be aliquoted and stored at -20 °C or lower.

Related Products:

Name	Cat #	Size
[Pyr3]- β -Amyloid (3-40)	AS-29906-01	0.1 mg
Pyr-FRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVV	AS-29906-1	1 mg
β -Amyloid (1-40)	AS-24235	0.5 mg
DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVV	AS-24236	1 mg
	AS-24236-5	5 mg

For Research Use Only