



## Product Information Sheet

---

Product Name:	Goat Anti-Rabbit IgG (H+L), highly cross-adsorbed, Hilyte Fluor™ 647-labeled
Catalog Number:	AS-61056-05-H647
Size:	0.5 mg
Concentration:	1mg/mL
Degree of Substitution: (DOS)	Dye/protein molar ratio is labeled on each vial.
Fluorescence:	Near-infrared fluorescence. Excitation/Emission wavelength= 649 nm/ 674 nm
Storage buffer:	10 mM phosphate, 150 mM NaCl, pH 7.2, with bovine serum albumin and 2 mM sodium azide as preservatives.
Storage:	The dye-protein conjugate is stable for 2~3 months at 4°C. For long-term storage, you can either divide the solution into aliquots and store at -20°C. Avoid multiple thaw-freeze cycles. Or you can add equal volume of glycerol (ACS grade or above) and store the solution at -20°C without aliquoting. The product will be stable for 1 year at -20°C.
Instructions:	The recommended concentration for most immunofluorescent staining is 0.1-10 µg/mL. If the non-specific binding background is high, you may decrease the concentration to minimize the background. You may also centrifuge the conjugate briefly and use the supernatant only for staining. This IgG is highly adsorbed with minimal cross-reaction to human, mouse, and rat serum protein

For *in vitro* research use only.