

Product Data Sheet

Product Name:	Orexin B, human	
Catalog Number:	AS-65573 (1 mg)	Lot Number: See label on vial
Sequence:	Arg-Ser-Gly-Pro-Pro-Gly-Leu-Gln-Gly-Arg-Leu-Gln-Arg-Leu-Leu-Gln-Ala- Ser-Gly-Asn-His-Ala-Ala-Gly-Ile-Leu-Thr-Met-NH2	
	RSGPPGLQGRLQRLLQASGNHAAGILTM-NH2	
Molecular Weight:	2898.9	
Peptide Purity:	Peak Area by HPLC ≥95%	
Storage:	This peptide is shipped at ambient temperature. Upon receipt, store lyophilized peptide at -20° C or lower. Reconstituted peptide can be aliquoted and stored at -20° C or lower.	
Description:	Orexins are synthesized by posterolateral n bind to orexin-1 (OxR1) and orexin 2 (OxR2 to both OxR1 and OxR2 while Orexin B can Orexins are involved in the regulation of aro homeostasis, and autonomic central control	e) receptors. Orexin A can bind associate with OxR2 only. Both busal and sleep states, energy
Related Products:	Orexin B, mouse, Cat# AS-65574 Orexin A, Cat# AS-24470	
References:	 Martynska L., et.al., <i>Neuroendocr Lett</i> 4 Ebrahim I.O., et.al., J of Royl Soc of Me Ramanjaneya M., et.al., <i>Endocrinology</i> 	d 95 , 227-230 (2002)

For Research Use Only