

IL2 ELISA Kit (Interleukin 2)

Antigen:	Interleukin 2 (IL2)
Reactivity:	Human
Application:	ELISA
Catalog no.:	ABIN1446208
Quantity:	5 x 96 tests

Product Details

Detection Method:	Colorimetric
Components:	Coat it Yourself ELISA. Capture Antibody, Biotinylated Detection Antibody, Standard, Streptavidin HRP and TMB.

Target details

Alternative Name:	IL-2 (IL2 ELISA Kit Abstract)
Research Area:	Cytokines, Immunology, Virology, Inflammation, Cancer
Pathways:	JAK-STAT Signaling

Application Details

Comment:	ELISA Sets including standard vials and essential components
Restrictions:	For Research Use only

Handling

Storage:	4 °C
Storage Comment:	Declared stable after 7 days at room temperature. Do not freeze.
Expiry Date:	1 week