

Datasheet for ABIN6952532

Coronavirus IgG Antibody ELISA Kit

Overview

Antigen: Coronavirus IgG Antibody

Reactivity: Rat Coronavirus (RCoV)

Method Type: Sandwich ELISA

Application: ELISA

Catalog no.

ABIN6952532

Coronavirus IgG Antibody ELISA Kit

Price

Please check our website for current price.

Quantity

2 x 96 tests

Availability

Available

Delivery in 3 to 4 Business Days (Ships to: Asia)

Product Details

Purpose: ELISA for the Qualitative determination of Anti-Coronavirus IgG in rat sera

Sample Type: Serum

Analytical Method: Qualitative

Detection Method: Colorimetric

Characteristics: Enzyme Immunoassay for the Qualitative determination of IgG Antibodies to Rat Coronavirus (SDA/RCV) in rat sera. Precoated plates using coronavirus protein for capture/detection

Target Details

Antigen: Coronavirus IgG Antibody

Application Details

Application Notes: Optimal working dilution should be determined by the investigator.

Plate: Pre-coated

Restrictions: For Research Use only

Handling

Storage: 4 °C

Please address orders & technical inquiries to:

antibodies-online GmbH
Schloss-Rahe-Str. 15
52072 Aachen
Germany

Phone +49 (0)241 95 163 153
Fax +49 (0)241 95 163 155
info@antibodies-online.com

If you prefer to order via fax you may use this as order form. Just enter all fields possible and send it to us using the fax number above.

Orderer

Salutation (Mr./Ms.):

Title:

Name:

Institution:

Phone:

Email:

Shipping address

Name recipient:

Institution:

Street:

Postal code:

City:

Country:

Billing address

Name recipient:

Institution:

Street:

Postal code:

City:

Country:

Other

VAT ID:

Purchase Order Number (PON) (Optional):

Your Order

Pos.	Product-ID	Product-Name	Amount
1.	ABIN6952532	Coronavirus IgG Antibody ELISA Kit	
2.			
3.			

I took note of and accepted the terms and conditions of antibodies-online GmbH (<https://www.antibodies-online.com/terms/>).

Place, Date

Signature