



Axygen® 200µL Universal Fit Pipet Tips, Bevelled, Graduated, Yellow, Nonsterile, 96 Tips/Rack, 5 Racks/Stack, 10 Stacks/Case (Product #TR-222-Y-STK)

Axygen® Universal Tips are designed to fit a wide variety of single and multi-channel pipettors. These polypropylene tips are nonpyrogenic and certified to be RNase-/DNase-free. They are available in Maxymum Recovery, bevelled, Eppendorf-style, nonsterile, sterile, Non-Filtered, and other features for specialized applications. Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject. Packaging includes, bulk pack, 96 tips/rack, and stacked rack.

- These polypropylene tips are nonpyrogenic and certified to be RNase-/DNase-free.
- Nonpyrogenic
- Graduation marks at 10, 50, and 100µL
- Bevelled



Product Details At-a-Glance

Qty./Pk.:	480
Qty./CS.:	4800
Capacity:	200 µL
Volume Range:	1-200 µL
Color:	Clear
Sterile:	No
Packaging Format:	96 Tips/Rack, 5 Racks/Stack, 10 Stacks/Case