

Anti-Mouse MHC Class II (I-A/I-E) PE

Catalogue Number : 86212-60

RUO: For Research Use Only. Not for use in diagnostic procedures.

Product Information

Clone: M5/114.15.2

Format/Conjugate: PE

Concentration: 0.2 mg/mL

Reactivity: Mouse

Laser: Blue (488nm), Yellow/Green (532-561nm)

Peak Emission: 578nm

Peak Excitation: 496nm

Filter: 585/40

Brightness (1=dim,5=brightest): 5

Isotype: Rat IgG2b, kappa

Formulation: Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH7.2.

Storage: Product should be kept at 2-8°C and protected from prolonged exposure to light.

Applications: FC

Description

The M5/114.15.2 monoclonal antibody or (M5/114) is specific for the mouse major histocompatibility complex class II, specific for the I-A and I-E subregion-encoded glycoproteins (alloantigens I-Ab, I-Ad, I-Aq, I-Ed, and I-Ek, but not the I-Af, I-Ak, nor I-As antigens). The antibody also reacts with cells from the H-2p or H-2r haplotypes mice.

The polymorphic determinant represented by the mouse major histocompatibility complex class II (MHC Class II) is expressed by the immune cells bearing the specified alloantigens, like T and B lymphocytes, monocytes, dendritic cells, and macrophages. The antibody doesn't react with the NOD H-2g7 haplotype alloantigens.

M5/114 mAb inhibits I-A-restricted T cell responses of the Haplotypes H-2b, H-2d, H-2q, and H-2u.

Preparation & Storage

The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified utilizing affinity chromatography and unreacted dye was removed from the product.

Application Notes

The antibody has been analyzed for quality through the flow cytometric analysis of the relevant cell type. For flow cytometric staining, the suggested use of this reagent is ≤0.06 ug per million cells in 100 µl volume. It is recommended that the reagent be titrated for optimal performance for each application.

References

1. Bhattacharya, A. L. O. K., Dorf, M. E., ; Springer, T. A. (1981). A shared alloantigenic determinant on Ia antigens encoded by the IA and IE subregions: evidence for I region gene duplication. *The Journal of Immunology*, 127(6), 2488-2495.
2. Mendiratta, S. K., Singh, N., Bal, V., ; Rath, S. (1996). Analysis of T-cell hybridomas with an unusual MHC class II-dependent ligand specificity. *Immunology*, 89(2), 238-244.
3. Unternaehrer, J. J., Chow, A., Pypaert, M., Inaba, K., ; Mellman, I. (2007). The tetraspanin CD9 mediates lateral association of MHC class II molecules on the dendritic cell surface. *Proceedings of the National Academy of Sciences*, 104(1), 234-239.