

bs-11181R-PE

[Conjugated Primary Antibody]



www.biossusa.com
support@biossusa.com
800.501.7654 [DOMESTIC]
+1.781.569.5821 [INTERNATIONAL]

Neurologin 3 Polyclonal Antibody, PE Conjugated

DATASHEET

Host: Rabbit

Target Protein: Neurologin 3

Clonality: Polyclonal

Isotype: IgG

Source: KLH conjugated synthetic peptide derived from human Neurologin 3

Purification: Purified by Protein A.

Storage: Aqueous buffered solution containing 100ug/ml BSA, 50% glycerol and 0.09% sodium azide.
Store at 4°C for 12 months.

Background: This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. Mutations in this gene may be associated with autism and Asperger syndrome. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Oct 2009].

Conjugation: PE

Excitation/ Emission: 480,565nm/578nm

Size: 100ul

Concentration: 1ug/ul

Applications: FCM(1:20-100)

Cross Reactive Species: Human
Mouse
Rat

For research use only. Not intended for diagnostic or therapeutic use.