

**bs-13547R-PE**

**[Conjugated Primary Antibody]**



www.biossusa.com  
support@biossusa.com  
800.501.7654 [DOMESTIC]  
+1.781.569.5821 [INTERNATIONAL]

ZAN Polyclonal Antibody, PE Conjugated

**DATASHEET**

**Host:** Rabbit

**Target Protein:** ZAN

**Clonality:** Polyclonal

**Isotype:** IgG

**Entrez Gene:** [7455](#)

**Source:** KLH conjugated synthetic peptide derived from human ZAN

**Purification:** Purified by Protein A.

**Storage:** Aqueous buffered solution containing 100ug/ml BSA, 50% glycerol and 0.09% sodium azide. Store at 4°C for 12 months.

**Background:** This gene encodes a sperm membrane protein that binds the zona pellucida of the egg in a species-specific manner. The encoded protein may be involved in signaling or gamete recognition. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008].

**Conjugation:** PE

**Excitation/ Emission:** 480,565nm/578nm

**Size:** 100ul

**Concentration:** 1ug/ul

**Applications:** FCM(1:20-100)

**Cross Reactive Species:** Human

For research use only. Not intended for diagnostic or therapeutic use.