

**bs-6282R-A647**

[Conjugated Primary Antibody]



www.biossusa.com  
support@biossusa.com  
800.501.7654 [DOMESTIC]  
+1.781.569.5821 [INTERNATIONAL]

Rabbit Anti-SERPINB4 Polyclonal Antibody, ALEXA FLUOR® 647  
Conjugated

— D A T A S H E E T —

**Host:** Rabbit

**Target Protein:** SERPINB4

**Clonality:** Polyclonal

**Isotype:** IgG

**Entrez Gene:** [23244](#)

**Swiss Prot:**

**Source:** KLH conjugated synthetic peptide derived from human SERPINB4

**Purification:** Purified by Protein A.

**Storage:** Aqueous buffered solution containing 100ug/ml BSA, 50% glycerol and 0.09% sodium azide.  
Store at 4°C for 12 months.

**Background:** May act as a protease inhibitor to modulate the host immune response against tumor cells.

**Conjugation:** ALEXA FLUOR® 647

**Excitation/ Emission:** 652nm/668nm

**Size:** 100ul

**Concentration:** 1ug/ul

**Applications:** IF(IHC-P)(1:50-200)

**Cross Reactive Species:** Human  
Mouse  
Rat

Caution: For research use only. Not for human or animal therapeutic or diagnostic use.