

bs-6282R-Cy5

[Conjugated Primary Antibody]

Bioss
ANTIBODIES

www.biossusa.com
support@biossusa.com
800.501.7654 [DOMESTIC]
+1.781.569.5821 [INTERNATIONAL]

Rabbit Anti-SERPINB4 Polyclonal Antibody, Cy5 Conjugated

— D A T — A S H — E E — T —

Host: Rabbit

Target Protein: SERPINB4

Clonality: Polyclonal

Isotype: IgG

Entrez Gene: [23244](#)

Swiss Prot:

Source: KLH conjugated synthetic peptide derived from human SERPINB4

Purification: Purified by Protein A.

Storage: Aqueous buffered solution containing 100ug/ml BSA, 50% glycerol and 0.09% sodium azide.
Store at 4°C for 12 months.

Background: May act as a protease inhibitor to modulate the host immune response against tumor cells.

Conjugation: Cy5

Excitation/ Emission: 625,650nm/670nm

Size: 100ul

Concentration: 1ug/ul

Applications: IF(IHC-P)(1:50-200)

Cross Reactive Species: Human
Mouse
Rat

Caution: For research use only. Not for human or animal therapeutic or diagnostic use.