

bs-6282R-HRP

[Conjugated Primary Antibody]



www.biossusa.com
support@biossusa.com
800.501.7654 [DOMESTIC]
+1.781.569.5821 [INTERNATIONAL]

Rabbit Anti-SERPINB4 Polyclonal Antibody, HRP Conjugated

— D A T — A — S — H — E — E — T —

Host: Rabbit

Target Protein: SERPINB4

Clonality: Polyclonal

Isotype: IgG

Entrez Gene: [23244](#)

Swiss Prot:

Source: KLH conjugated synthetic peptide derived from human SERPINB4

Purification: Purified by Protein A.

Storage: Aqueous buffered solution containing 100ug/ml BSA, 50% glycerol and 0.09% Gentamicin.
Store at 4°C for 12 months.

Background: May act as a protease inhibitor to modulate the host immune response against tumor cells.

Conjugation: HRP

Excitation/ Emission: N/A

Size: 100ul

Concentration: 1ug/ul

Applications: WB(1:100-1000)
IHC-P(1:100-500)

Cross Reactive Species: Human
Mouse
Rat

Caution: For research use only. Not for human or animal therapeutic or diagnostic use.