

**bs-6282R**

[Primary Antibody]



www.biossusa.com  
support@biossusa.com  
800.501.7654 [DOMESTIC]  
+1.781.569.5821 [INTERNATIONAL]

## SERPINB4 Polyclonal Antibody

### DATASHEET

**Host:** Rabbit

**Target Protein:** SERPINB4

**Clonality:** Polyclonal

**Isotype:** IgG

**Entrez Gene:** [23244](#)

**Source:** KLH conjugated synthetic peptide derived from human SERPINB4

**Purification:** Purified by Protein A.

**Storage:** Aqueous buffered solution containing 100ug/ml BSA, 50% glycerol and 0.09% sodium azide.  
Store at -20 °C for 12 months.

**Background:** May act as a protease inhibitor to modulate the host immune response against tumor cells.

**Size:** 100ul

**Concentration:** 1ug/ul

**Applications:** WB(1:100-1000)  
IHC-P(1:100-500)  
IF(IHC-P)(1:50-200)

**Cross Reactive:** Human

**Species:** Mouse  
Rat

For research use only. Not intended for diagnostic or therapeutic use.

### VALIDATION IMAGES