



80 Webster Street  
Worcester, MA 01603 (USA)  
TEL: 508-770-1190

Web Site: BioPAL.com  
e-mail: cs@BioPAL.com

## Product Data Sheet

---

Product Name:	Molday ION C6Amine
Catalog Number:	CL-50Q02-6A
Description:	Stable aqueous suspension of iron-based superparamagnetic nanoparticles (USPIO)
Application:	Used to label cells, providing MRI contrast, magnetic targeting, capture and retention.
Packaging:	2 ml sealed serum bottle
Volume:	2.0 ml
Weight of Iron:	5 mg Fe/ml
Mean Particle Size:	35 nm
Core:	Magnetite
Mean Core Size:	8 nm
Coating:	Dextran
Number of Particles:	~6.4E16 per g
Density:	~1.25 g/cm <sup>3</sup>
Zeta Potential:	~ +48 mV
R1 (mM <sup>-1</sup> sec <sup>-1</sup> )	32.1
R2 (mM <sup>-1</sup> sec <sup>-1</sup> )	67.1
Cell Labeling:	Universal cell labeling reagent - No transfection agent required
Functionality	Amine
Visualization:	Prussian blue, MRI, TEM
Storage Buffer:	dH <sub>2</sub> O
Autoclaved:	No
Long-Term Storage:	4 °C - <b>Do not freeze!</b>
Similar Products:	None
Expiration Date:	None

---

### Known References:

---

Curr Gene Ther. 2014;14(2):136-45.  
Int J Nanomedicine. 2014 Jan 8;9:337-50  
PLoS One. 2013 Sep 13;8(9):e74658  
Mol Med Rep. 2013 Nov;8(5):1446-52  
Radiology. 2012 Oct;265(1):175-85  
PLoS One. 2011;6(9):e24730  
Magn Reson Med. 2011 Feb;65(2):564-74  
Biotechniques. 2011 Apr;50(4):223-7  
Contrast Media Mol Imaging. 2011 Nov-Dec;6(6):514-22  
Contrast Media Mol Imaging. 2011 Jan-Feb;6(1):7-18  
Sci China Life Sci. 2011 Nov;54(11):981-7  
J Am Coll Cardiol. 2009 Oct 20;54(17):1627-8