



Mouse anti-c-Myc IgG conjugated to Alkaline Phosphatase

Product Number AP-1714
Amount 100 μ g
Store at 4°C

Form/Shipping & Storage

Supplied as a lyophilized powder. Upon receipt, store at 4°C in the dark. Phycobiliproteins are sensitive to freeze-thaw cycles; after reconstitution, store at 2-8°C in the dark, do not freeze.

Handling

Prior to use, reconstitute to 1 ml with 50% glycerol, vortex gently and allow material to reconstitute for 20 minutes. Avoid exposure to heat.

Buffer

Upon rehydration with 1 ml of 50% glycerol, the product is in 25mM Tris (pH 7.2), 150mM sodium chloride, and 2mM sodium azide as a preservative.

Stability

Lyophilized material is stable for up to one year. After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 6 months. If further dilution of the conjugate is required, use diluted material within one week.

Antibody/Antigen Information

Mouse IgG₁ clone 9E10 against sequence EQKLISEEDL.

Note

For research use only, not for diagnostic or therapeutic use.



COLUMBIA BIOSCIENCES