



# PRODUCT Hu-A006 <u>ANTI-HUMAN CYP4F+ IgG</u> Polyclonal Antibody Developed in Rabbits, IgG Fraction LOT RaP/B#4+5

Antiserum was developed in rabbits using purified human liver CYP4F2 as immunogen. The whole IgG fraction was purified from antiserum using caprylic acid/ammonium sulfate fractionation. Anti-human CYP4F+ IgG is provided in a lyophilized state after freeze-drying from 100 mM potassium phosphate buffer (pH 7.4), 150 mM KCI, and 2.5 µM thimerosal (added as a preservative).

# • Specificity and Purity

Specificity has been determined on Western blots (see below). Anti-human CYP4F+ IgG reacts strongly with CYP4F2 and several other 56-57 kDa proteins found in human liver microsomes. <u>Due to extensive structural homology among CYP4F2 and other human CYP4F P450s</u>, significant cross-reactivity of anti-CYP4F+ IgG with CYP4F3a, CYP4F3b, CYP4F11 and CYP4F12 is observed. A cross-reaction with a 70 kDa non-P450 protein is also sometimes noted, depending upon the liver sample being tested. Antibody purity has been established by denaturing SDS-PAGE, giving two protein bands with  $M_R$  of 50 kDa and 25 kDa as visualized Coomassie blue staining, which correspond to the heavy and light chains, respectively, of rabbit IgG.

# <u>Reconstitution of Lyophylized Product and Storage</u>

Store the lyophylized product at 0-5°C. For Western blotting, reconstitute by adding 1 ml of PBS/50% glycerol to one vial of lyophilized IgG (1 mg) and mix vial gently until powder dissolves. Upon reconstitution, the solution can be stored at -20°C, as the presence of glycerol will prevent freeze/thaw cycles. Anti-CYP4F+ IgG solutions without glycerol should be also be stored at -20°C but subjected to freeze/thaw cycles as seldom as possible.



# Immunochemical Characteristics of Anti-CYP4F+ IgG

Lane A = Liver Microsomes from Subject K (15  $\mu$ g) Lane B = rCYP4F2 expressed in *Sf9* insect cell lysates (0.5  $\mu$ g) Lane C = rCYP4F3b expressed in *Sf9* insect cell lysates (0.5  $\mu$ g) Lane D = rCYP4F11 expressed in *Sf9* insect cell lysates (0.5  $\mu$ g) Lane E = rCYP4F12 expressed in *Sf9* insect cell lysates (1.0  $\mu$ g) Lane F = rCYP4F3a expressed in *T. ni* insect cells (0.5  $\mu$ g) Lane G = Liver Microsomes from Subject L (15  $\mu$ g)

### Use for Western Blotting

Incubate blots overnight with 2.5-5.0 µg rabbit anti-human CYP4F+ IgG IgG/ml of appropriate blocking solution. After washing to remove unbound antibody, incubate with anti-rabbit IgG conjugate of choice (e.g, anti-rabbit IgG-peroxidase or anti-rabbit IgG-biotin), and develop accordingly. A detailed Western blotting method can be found in the <u>PROTOCOLS</u> section.

#### • Use for Immunoinhibition

Incubation of anti-human CYP4F+ IgG with human liver microsomes at a ratio of 1.7 mg IgG/mg microsomal protein before reaction initiation will typically give 90% inhibition of exemplary CYP4F-catalyzed reactions [e.g., arachidonate or leukotriene B4 (LTB4) ω-hydroxylation; see below]. The methodology for conducting P450 immunoinhibition assays is given in the <u>PROTOCOLS</u> section.

