
Gold nanoparticles - 60nm

Product details

Product Code: EM.GC60

Type: Gold Nanoparticles

Applications: Life Science Lateral Flow

Product summary

Optical Density: 1.19 @ 520nm

Presentation matrix: Suspended in H₂O, no preservative Residual chemical left from manufacture

Specification: 60 nm

Application notes

Recommended Use: In vitro research

Tested Applications: Lateral flow

Additional information

Storage: Store at 2-8°C – do not freeze

Shelf Life: 12 months

Keywords: gold colloid, gold nanoparticles
