



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06009 - Goat Anti-FOXD2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FOXD2 antibody, forkhead box D2 antibody, FKHL17 antibody, FREAC9 antibody, freac-9 antibody, forkhead-related activator 9 antibody, forkhead (Drosophila)-like 17 antibody, Forkhead, drosophila, homolog-like 17 antibody, FREAC-9 antibody, forkhead-like 17 antibody

Official Symbol: FOXD2

Accession Number(s): NP_004465.3

Human GeneID(s): [2306](#)

Non-Human GeneID(s): 17301 (mouse)

Immunogen

Peptide with sequence FASKVAGLSGCHF, from the C Terminus of the protein sequence according to NP_004465.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: No signal obtained yet but low background observed in Human Kidney and 293 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat