



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07582 - Goat Anti-AOC3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: AOC3 antibody, amine oxidase, copper containing 3 (vascular adhesion protein 1) antibody, HPAO antibody, SSAO antibody, VAP-1 antibody, VAP1 antibody, amine oxidase, copper containing 3 antibody, copper amine oxidase antibody, semicarbazide-sensitive amine oxidase antibody, vascular adhesion protein 1 antibody

Official Symbol: AOC3

Accession Number(s): NP_003725.1

Human GeneID(s): [8639](#)

Non-Human GeneID(s): 11754 (mouse), 29473 (rat)

Immunogen

Peptide with sequence C-DEDPSFYSDSIY, from the internal region (near the C Terminus) of the protein sequence according to NP_003725.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.


Western blot: Preliminary experiments gave an approx 55kDa band in Human Adipose and Lung lysates after 0.01µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 84.6kDa according to NP_003725.1. The 55kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

IHC: In paraffin embedded Human Liver shows membranous staining in sinusoid lining cells surrounding hepatocytes. Recommended concentration, 3-5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat



EB07582 (3.8µg/ml) staining of paraffin embedded Human Liver. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.