
Anti-Medaka vitellogenin IgG,
Horseradish Peroxidase conjugate

Code : MVD33-HRP

Lot :

Lot size : 0.5 mL

Background: Vitellogenin, the egg yolk precursor protein, is induced by estrogens.

On the other hand, environmental estrogens are known to disrupt endocrine system in animals. Recently, vitellogenin has proven to be an ideal marker for enviromental estrogens.

Immunogen Purified medaka vitellogenin

Source :

Applications : ELISA, Western blot

Normal dilution ELISA : 1:100 ~ 1:1000

range : Western blot : 1:100 ~ 1:5000

Storage : Short term, 4°C

Long term, -20°C

Avoid freeze thaw cycles

Note : Since assay conditions vary, the optimum dilution should be determined for particular application.

Warning : For research use only

FUJIKURA KASEI CO.,LTD.

R&D Center

13-1, Sakurada, 5-chome, Kuki-shi, Saitama, 340-0203 Japan

Tel: +81-480-57-1155 E-mail: ENBIO@fkkasei.co.jp