




Myc (Ab-373) Antibody

E021036

- Catalog Number:** E021036-1, E021036-2
Amount: 50µg/50µl, 100µg/100µl
Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage/Stability: Store at -20°C/1 year
Immunogen: The antiserum was produced against synthesized non-phosphopeptide derived from human Myc around the phosphorylation site of serine 373 (K-R-S^P-F-F).
Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Specificity/Sensitivity: Myc (Ab-373) antibody detects endogenous levels of total Myc protein.
Reactivity: Human, Mouse, Rat
Applications: WB: 1:500~1:1000
Swiss-Prot No. : P01106
References: Baudino T A, et al. (2001) Mol Cell Biol. 21: 691-702.
Blackwood E M, et al. (1991) Science. 251:1211-1217.
Henriksson M, et al. (1996) Adv Cancer Res. 68: 109-182.
Grandori C, et al. (2000) Annu Rev Cell Dev Biol. 16: 653-699.



For Research Use Only