




SAPK/JNK (Ab-183) Antibody


E021241

- Catalog Number:** E021241-1, E021241-2
- Amount:** 50µg/50µl, 100µg/100µl
- Swiss-Prot No. :** P45984; P53779
- Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
- Storage/Stability:** Store at -20°C/1 year
- Immunogen:** The antiserum was produced against synthesized non-phosphopeptide derived from human SAPK/JNK around the phosphorylation site of threonine 183 (M-M-T^P-P-Y).
- Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
- Specificity/Sensitivity:** SAPK/JNK (Ab-183) Antibody detects endogenous levels of total SAPK/JNK protein.
- Reactivity:** Human, Mouse, Rat
- Applications:** WB: 1:500~1:1000 IHC: 1:50~1:100
- References:** I. Ferrer, et al. (2003) *Neuropathology & Applied Neurobiology* 29: 23
 Zhonghong Guan, et al. (1999) *J Biol Chem*, Vol. 274: 36200-36206
 D.Margriet Ouwens1, et al. (2002) *The EMBO Journal* 21: 3782-3793,
 S Matsuyoshi1, et al. (2006) *British Journal of Cancer* 94: 532



Peptide - +

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using SAPK/JNK (Ab-183) Antibody (E021241).



Western blot analysis of extracts from 293 cell using SAPK/JNK (Ab-183) Antibody (E021241, Lane 1, 2) and SAPK/JNK (phospho-Thr183) antibody (E011249, Lane3, 4).

For Research Use Only