




OCT3 Antibody

E11-0283C

- Catalog Number:** E11-0283C
- Amount:** 100µg/100µl
- Swiss-Prot No. :** Q01860
- All Names:** OCT3, OCT4, OTF3, OTF4, OCTAMER-BINDING TRANSCRIPTION FACTOR 3, OCTAMER-BINDING TRANSCRIPTION FACTOR 4
- Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
- Storage/Stability:** Store at -20°C/1 year
- Immunogen:** The antiserum was produced against synthesized peptide derived from human OCT3 antibody.
- Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
- Specificity/Sensitivity** OCT3 antibody detects endogenous levels of total OCT3 protein.
- Reactivity:** Human, Mouse, Rat
- Applications:** WB: 1:500~1:1000 IHC: 1:50~1:100
ELISA: 1:10000
- References:** Y Saijoh, Genes Cells, Feb 1996; 1: 239 - 252.
Christopher L. Shaffer, Drug Metab. Dispos., Apr 2007; 10:1124.
Brian Boer, Nucleic Acids Res., Mar 2007; 35: 1773 - 1786.
E Pikarsky, Mol. Cell. Biol., Feb 1994; 14: 1026 - 1038.



For Research Use Only

For Research Use Only