



N-Myc Antibody

E300154

Product description:

Formulation: Lyophilized anti-N-Myc

Antibody serum containing 0.02% NaN_3

Catalog No.: E300154

Size: 100ug/200ul

Source:

Rabbit polyclonal N-Myc⁽¹⁾ antibody was raised against a 14 amino acid peptide at the intermediate residues of human N-Myc protein.

Species Reactivity:

Human, mus

Not yet tested in other species.


Application:

For research use only

Peptide ELISA titer: 1:80,000

Western blot analysis: 1:500-1,000

Optimal dilutions/concentrations should be determined by the end user.



Dilution instruction and storage:

Already reconstituted in 200ul sterile H_2O . The antibody is stable for at least one year when stored at -20°C after reconstitution.

Avoid repeated freeze/thaw cycles.

Reference:

(1) Prendergast, G.C., Lawe, D., and Ziff, E.B. 1991. Association of Myc, the murine homolog of Max, with c-Myc stimulates methylation-sensitive DNA binding and Ras cotransformation. Cell.1991; 65(3):395-407.

For Research Use Only