


EZRIN

Gene-tagged cell line (HeLa)

Catalog no: EXP-001



Cell type: HeLa

Gene symbol: EZR

NCBI gene ID: 7430

Protein: Ezrin

Subcellular location: Cell membrane

Modification: C-terminal mRuby3

Excitation / Emission (nm): 558 / 592

Antibiotic resistance: ZeocinTM

Population type: Homozygous

Protein summary from NCBI database


As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene. [provided by RefSeq, Jul 2008]

info@xpresscells.com www.xpresscells.com

CUSTOM CELL LINE SERVICES AVAILABLE UP TO 3 KNOCK-INS IN A SINGLE CELL LINE

ExpressCells' FAST-HDR knock-in technology

ExpressCells uses CRISPR and FAST-HDR vector technology to knock-in fluorescent, luminescent, or other tags at the C-terminus of endogenous genes. The non-viral FAST-HDR system enables rapid labeling of up to three proteins of interest in a single mammalian cell line.



Handling

Culture medium: Dulbecco's Modified Eagle Medium (DMEM), high glucose supplemented with 10% fetal bovine serum (FBS) and penicillin/streptomycin to prevent bacterial contamination.

Thawing: Transfer the frozen tube to a 37° C water bath and let contents thaw. Transfer tube contents to 10 mL of prewarmed medium in a biosafety hood and centrifuge at $200 \times g$ for 5 min. Resuspend the pellet in 5 mL of medium and transfer to a mammalian cell culture flask.

Safety: Biosafety level 2.

References

- 1. Gene [database online]. Washington DC: NCBI; 2020. https://www.ncbi.nlm.nih.gov/gene/7430. Accessed February 6, 2020.
- 2. Bretsche A. Regulation of cortical structure by the ezrin-radixin-moesin protein family. *Curr Opin Cell Biol.* 1999;11:109-116.
- Jiang Wang, Xiaolin Zhao, Jun Qi, Caihong Yang, Hao Cheng, Ye Ren, Lei Huan. Eight proteins play critical roles in RCC with bone metastasis via mitochondrial dysfunction. Clin Exp Metastasis. 2015;32:605-622.
- Perez-Leal O, Nixon-Abell J, Barrero CA, Gordon J, Rico MC. A versatile vector system for the fast generation of knock-in cell lines with CRISPR [preprint published online February 6 2020]. bioRxiv. doi: 10.1101/2020.02.06.927384.

For research use only.