

11-712-C025

## Monoclonal Antibody to CD3 zeta (Phospho-Tyr72) Purified Antibody (0.025 mg)

<b>Clone:</b>	EM-26
<b>Isotype:</b>	Mouse IgG2b
<b>Specificity:</b>	The mouse monoclonal antibody EM-26 recognizes phosphorylated tyrosine 72 of CD3 zeta chain (CD247), which is a component of TCR/CD3 complex expressed on T cells.
<b>Regulatory Status:</b>	RUO
<b>Immunogen:</b>	A phospho specific peptide corresponding to the amino acids surrounding tyrosine 72 of mouse CD3 zeta linked to KLH
<b>Species Reactivity:</b>	Human, Mouse
<b>Application:</b>	Flow Cytometry Recommended dilution: 2 - 10 µg/ml Positive control: Jurkat cells treated with pervanadate; T-cells from lymph nodes of OT-1 mouse treated with pervanadate Western Blotting Recommended dilution: 2 - 5 µg/ml Positive control: Jurkat cells lysate treated with pervanadate; Splenocyte lysate of Balb/c or F1 mouse treated with pervanadate Application note: Non-reducing conditions recommended Immunocytochemistry
<b>Purity:</b>	> 95% (by SDS-PAGE)
<b>Purification:</b>	Purified by protein-A affinity chromatography
<b>Concentration:</b>	1 mg/ml
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
<b>Storage / Stability:</b>	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.
<b>Expiration:</b>	See vial label
<b>Lot Number:</b>	See vial label
<b>Background:</b>	CD3 complex is crucial in transducing antigen-recognition signals into the cytoplasm of T cells and in regulating the cell surface expression of the TCR complex. T cell activation through the antigen receptor (TCR) involves the cytoplasmic tails of the CD3 subunits CD3 gamma, CD3 delta, CD3 epsilon and CD3 zeta (CD247). These CD3 subunits are structurally related members of the immunoglobulins super family encoded by closely linked genes on human chromosome 11. The CD3 components have long cytoplasmic tails that associate with cytoplasmic signal transduction molecules. This association is mediated at least in part by a double tyrosine-based motif present in a single copy in the CD3 subunits. CD3 may play a role in TCR-induced growth arrest, cell survival and proliferation.

**For laboratory research only, not for drug, diagnostic or other use.**



**Antibodies**

- References:**
- \*Dopfer EP, Schöpf B, Louis-Dit-Sully C, Dengler E, Höhne K, Klescová A, Prouza M, Suchanek M, Reth M, Schamel WW: Analysis of novel phospho-ITAM specific antibodies in a S2 reconstitution system for TCR-CD3 signalling. Immunol Lett. 2010 May 4;130(1-2):43-50
  - \*Borroto A, Arellano I, Dopfer EP, Prouza M, Suchanek M, Fuentes M, Orfao A, Schamel WW, Alarcón B: Nck recruitment to the TCR required for ZAP70 activation during thymic development. J Immunol. 2013 Feb 1;190(3):1103-12.

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic or therapeutic use. Not for resale or transfer either as a stand-alone product or as a component of another product without written consent of EXBIO. EXBIO will not be held responsible for patent infringement or other violations that may occur with the use of our products. All orders are accepted subject to EXBIO's term and conditions which are available at [www.exbio.cz](http://www.exbio.cz).

**For laboratory research only, not for drug, diagnostic or other use.**

---

EXBIO Praha | Nad Safinou II 341 | 252 50 Vestec u Prahy | Czech Republic  
Tel: +420 261 090 666 | Fax: +420 261 090 660 | [orders@exbio.cz](mailto:orders@exbio.cz) | [www.exbio.cz](http://www.exbio.cz)