

1F-626-C100

Monoclonal Antibody to CD62L / L-Selectin (mouse) Fluorescein (FITC) conjugated (0.1 mg)

Clone:	MEL-14
Isotype:	Rat IgG2a
Specificity:	The rat monoclonal antibody MEL-14 reacts with mouse CD62L (L-selectin), a 75 kDa single chain type I glycoprotein expressed on most peripheral blood B lymphocytes, T lymphocytes, monocytes and granulocytes; it is also present on a subset of NK cells and certain hematopoietic malignant cells.
Regulatory Status:	RUO
Immunogen:	C3H/eb mouse B cell lymphoma 38C-13
Species Reactivity:	Mouse
Preparation:	The purified antibody is conjugated with Fluorescein isothiocyanate (FITC) under optimum conditions. The reagent is free of unconjugated FITC.
Concentration:	0.5 mg/ml
Storage Buffer:	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage / Stability:	Store in the dark at 2-8°C. Do not freeze. Avoid prolonged exposure to light. Do not use after expiration date stamped on vial label.
Usage:	The reagent is designed for Flow Cytometry analysis. Suggested working concentration is 6 µg/ml. Indicated dilution is recommended starting point for use of this product. Working concentrations should be determined by the investigator.
Expiration:	See vial label
Lot Number:	See vial label
Background:	CD62L (L-selectin) is an adhesion glycoprotein that is constitutively expressed on the cell surface of leukocytes and mediates their homing to inflammatory sites and peripheral lymph nodes by enabling rolling along the venular wall. CD62L is also involved in activation-induced neutrophil aggregation. Activation-dependent CD62L shedding, however, counteracts neutrophil rolling. CD62L has also signaling roles including enhance of chemokine receptor expression. Similarly to CD62P, the major ligand of CD62L is PSGL-1 (P-selectin glycoprotein ligand-1).

For laboratory research only, not for drug, diagnostic or other use.

**Antibodies****References:**

- *Gallatin WM, Weissman IL, Butcher EC: A cell-surface molecule involved in organ-specific homing of lymphocytes. *Nature* 1983 Jul 7;304:30-34.
- *Brawand P, Cerottini JC, MacDonald HR: Hierarchical utilization of different T-cell receptor Vbeta gene segments in the CD8(+)-T-cell response to an immunodominant Moloney leukemia virus-encoded epitope in vivo. *J Virol.* 1999 Nov;73(11):9161-9.
- *Fisson S, Darrasse-Je G, Litvinova E, Septier F, Klatzmann D, Liblau R, Salomon BL: Continuous activation of autoreactive CD4+ CD25+ regulatory T cells in the steady state. *J Exp Med.* 2003 Sep 1;198(5):737-46.
- *ten Bruggencate SJ, Hillyer LM, Woodward BD: The proportion of CD45RA(+)/CD62L(+) (quiescent-phenotype) T cells within the CD8(+) subset increases in advanced weight loss in the protein- or energy-deficient weanling mouse. *J Nutr.* 2001 Dec;131(12):3266-9.
- *Mannering SI, Cheers C: Interleukin-2 and loss of immunity in experimental *Mycobacterium avium* infection. *Infect Immun.* 2002 Jan;70(1):27-35.
- *Pop SM, Wong CP, He Q, Wang Y, Wallet MA, Goudy KS, Tisch R: The type and frequency of immunoregulatory CD4+ T-cells govern the efficacy of antigen-specific immunotherapy in nonobese diabetic mice. *Diabetes.* 2007 May;56(5):1395-402.
- *Chin CS, Miller CH, Graham L, Parviz M, Zacur S, Patel B, Duong A, Bear HD: Bryostatin 1/ionomycin (B/I) ex vivo stimulation preferentially activates L-selectin-low tumor-sensitized lymphocytes. *Int Immunol.* 2004 Sep;16(9):1283-94.
- *Chen BJ, Cui X, Sempowski GD, Liu C, Chao NJ: Transfer of allogeneic CD62L-memory T cells without graft-versus-host disease. *Blood.* 2004 Feb 15;103(4):1534-41.
- *Xu H, Manivannan A, Crane I, Dawson R, Liversidge J: Critical but divergent roles for CD62L and CD44 in directing blood monocyte trafficking in vivo during inflammation. *Blood.* 2008 Aug 15;112(4):1166-74.
- *Friedline RH, Wong CP, Steeber DA, Tedder TF, Tisch R: L-selectin is not required for T cell-mediated autoimmune diabetes. *J Immunol.* 2002 Mar 15;168(6):2659-66.
- *Richards H, Longhi MP, Wright K, Gallimore A, Ager A: CD62L (L-selectin) down-regulation does not affect memory T cell distribution but failure to shed compromises anti-viral immunity. *J Immunol.* 2008 Jan 1;180(1):198-206.
- *And many other.

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic or therapeutic use. Not for resale or transfer either as a stand-alone product or as a component of another product without written consent of EXBIO. EXBIO will not be held responsible for patent infringement or other violations that may occur with the use of our products. All orders are accepted subject to EXBIO's term and conditions which are available at www.exbio.cz.

For laboratory research only, not for drug, diagnostic or other use.

EXBIO Praha | Nad Safinou II 341 | 252 50 Vestec u Prahy | Czech Republic
Tel: +420 261 090 666 | Fax: +420 261 090 660 | orders@exbio.cz | www.exbio.cz