

Combo: Clophosome®-A and Control Liposomes (anionic)



SKU: F70101C-AC

* Volume: 2 mL
 10mL



Description

Similar Products

Product Description

Combo: Clophosome®-A and Control Liposomes

Included in this convenience combo kit are 1 vial of the anionic clodronate liposomes and 1 vial of the anionic control liposomes. Your save is more than 20%.

Product code: F70101C-AC

This kit consists of the following items

1 vial (2 or 10mL) of Clophosome®-A - anionic clodronate liposomes ([F70101C-A](#))

1 vial (2 or 10mL) of control anionic liposomes ([F70101-A](#))

Store in refrigerator (2-8 degree C)