

MAB\_0010  
Product Data Sheet  
FOR RESEARCH PURPOSES ONLY . NOT FOR USE IN HUMANS

**Name:** Anti-ceramide  
**Clone:** S58-9  
**Immunogen:** not declared  
**Specificity:** ceramide  
**Storage:** room temperature (up to 7 d); -16 to -24°C (long term storage)  
**Type:** monoclonal  
**Species:** mouse  
**Isotype:** IgM  
**Preparation:** ascites  
**Form:** fluid  
**Amount:** 100  $\mu$ l, 200  $\mu$ l, 500  $\mu$ l  
**Concentration:** at least 1 mg/ml  
**Additives:** 0.02 % NaN<sub>3</sub>  
**Active in:** EIA, blot, IEM, IHC, IF  
**Description:** No cross-reactivity with cerebrosides, sphingomyelin, cholesterol or any other lipid when appropriate blocking conditions are choosen. For details see references.