

MAB_0011
Product Data Sheet
FOR RESEARCH PURPOSES ONLY . NOT FOR USE IN HUMANS

Name: Anti-ceramide
Clone: S58-9
Immunogen: not declared
Specificity: ceramide
Storage: room temperature (up to 7 d); -16 to -24°C (long term storage)
Type: monoclonal
Species: mouse
Isotype: IgM
Preparation: serum-free supernatant
Form: fluid
Amount: 200 μ l, 400 μ l, 1 ml
Concentration: at least 25 μ g/ml
Additives: 0.02 % NaN₃
Active in: EIA, blot, IEM, IHC, IF
Description: No cross-reactivity with cerebrosides, sphingomyelin, cholesterol or any other lipid when appropriate blocking conditions are choosen. For details see references.