

MAB_0013
Product Data Sheet
FOR RESEARCH PURPOSES ONLY . NOT FOR USE IN HUMANS

Name: Anti-ceramide
Clone: S58-9
Immunogen: not declared
Specificity: ceramide
Storage: room temperature (up to 7 d); -16 to -24°C (long term storage)
Type: monoclonal
Species: mouse
Isotype: IgM
Preparation: FCS miniperm supernatant
Form: fluid
Amount: 100 μ l, 200 μ l, 500 μ l
Concentration: at least 250 μ g/ml
Additives: 0.02 % NaN₃
Active in: EIA, blot, IEM, IHC, IF
Description: No cross-reactivity with cerebrosides, sphingomyelin, cholesterol or any other lipid when appropriate blocking conditions are choosen. For details see references.