

Human Papillomavirus type 18 E7 protein, His tag

Cat. No. GTX133412-pro**Species** Human papillomavirus type 18**Package**
100 µg


PROPERTIES

Form Liquid**Buffer** PBS, 0.1% SDS**Preservative** No preservative**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 3 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein of Human Papilloma virus type 18 E7 protein**Expression System** E. Coli**Purity** >95%**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX133412-pro Image**

GTX133412-pro Human Papillomavirus type 18 E7 protein was analyzed using SDS-PAGE and stained with coomassie blue and captured by monochrome camera.

For full product information, images and publications, please visit our [website](#).