



Product Data Sheet

| | |
|--|---------------------------|
| Product: | Human Factor IXa - EGR |
| Catalog number: | HCIXA-EGR |
| Concentration: | Lot dependent |
| Specific activity ¹ : | Functionally inactive |
| Extinction Coefficient: $E_{1\text{cm},280\text{ nm}}^{1\%} =$ | 14.0 |
| Molecular weight: | 45,000 |
| Standard buffer: | 50% Glycerol/Water (v/v)* |
| Recommended Storage: | -20°C |

PRODUCT NOTES:

1. Determined by assay for competitive inhibition with purified human IXa.

* This product may contain residual ammonium sulfate, and may require dialysis prior to use in applications that are sensitive to the presence of free amines.

CAUTION: Product is of human blood/plasma origin. Although the starting material has been tested and found negative or non-reactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material.

This product is to be used for **IN VITRO** research purposes only, and is not intended for clinical or diagnostic use!