



Product Data Sheet

Product:	Sheep anti-Human Factor VIII - Polyclonal
Catalog number:	PAHFVIII-S
Concentration:	Lot dependent
Specificity:	Factor VIII
Extinction Coefficient:	$E_{1\text{cm},280\text{ nm}}^{1\%} = 14.0$
Subclass:	IgG
Molecular weight:	150,000
Standard buffer:	50% Glycerol/HBS, pH 8.0 (v/v)*
Recommended Storage:	-20°C

PRODUCT NOTES:

Properties: Western blot (), ELISA (+), RIEP (); useful for purification

* This product may contain residual ammonium sulfate, and may require dialysis prior to use in applications that are sensitive to the presence of free amines.

CAUTION:

This product is to be used for **IN VITRO** research purposes only, and is not intended for clinical or diagnostic use!