



**CATALOGUE #:** 3E1

**PRODUCT NAME:** Monoclonal mouse anti-Ebola virus, Zaire strain

---

**MAbs:** FE18, FE25, FE37

**MAb isotypes:** IgG2a for all MAbs

**Applications:** All MAbs are working in ELISA.

MAbs FE18 and FE37 are working in Western blotting. MAbs FE18 and FE37 recognize high molecular weight material under non-reducing conditions and ~30K protein band under reducing conditions in Western blotting.

**Purification:** Chromatography on protein G Sepharose

**Presentation:** PBS, pH 7.2, 0.1 % sodium azide (NaN<sub>3</sub>)

**Storage:** + 4 °C

**Material safety note:** This product is sold as an antibody preparation **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Contains sodium azide (0.1 %) as preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.