



**CATALOGUE #:** 8BFP

**PRODUCT NAME:** BNP and NT-proBNP free plasma

---

- Source:** Pooled normal human plasma.  
Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies and HCV.
- Application:** BNP and NT-proBNP free plasma is free of whole proBNP molecule and any proBNP fragments (BNP and NTproBNP). It can be used for preparation of any proBNP-related standards and calibrators for immunoassays based on detection of any kinds of proBNP-related peptides.
- Purification:** Immunoaffinity chromatography of pooled normal human plasma.  
Gel for affinity column was prepared by coupling anti-BNP and anti-NTproBNP monoclonal antibodies with different epitope specificity with BrCN-activated Sepharose CL 4B.
- Presentation:** Frozen liquid.
- Storage:** - 20 °C