



CATALOGUE #: 8G16-17b

PRODUCT NAME: Disialoganglioside GD1a-NAcGal, bovine

Source: Bovine brain tissue.

Applications: One of the minor neuronal gangliosides. Has the terminal GalNAc-residue with beta linkage to the terminal Gal-residue of GD1a. It is stable against the majority of sialidases, including the sialidase from *Vibrio cholerae*.

Analysis: Purity more than 98 % by HPTLC. MW 2030

Presentation: Lyophilized from water solution.

Storage: -20 °C in lyophilized form.
After reconstitution the solution should be stored at -20 °C.