



CATALOGUE #: 8G16-2b

PRODUCT NAME: Monosialoganglioside GM1, bovine

Source: Bovine brain tissue.

Applications: Useful as an antigen and receptor for cholera toxin, growth inhibitor marker in fibroblasts and marker for lymphoid subpopulations. Stable for the majority of sialidase.

Analysis: Purity more than 98 % by HPTLC. MW 1545

Presentation: Lyophilized from water solution.

Storage: -20 °C in lyophilized form.
After reconstitution the solution should be stored at -20 °C.