



CATALOGUE #: 8G16-6b

PRODUCT NAME: Disialoganglioside GD1a, bovine

Source: Bovine brain tissue.

Applications: One of the main neuronal gangliosides. Useful as a differentiation marker of tumor cell growth. Can be converted to GM1 by bacterial and other sialidase from mammalian sources without detergent.

Analysis: Purity more than 98 % by HPTLC. MW 1827

Presentation: Lyophilized from water solution.

Storage: -20 °C in lyophilized form.
After reconstitution the solution should be stored at -20 °C.