



CATALOGUE #: 8IN74-1

PRODUCT NAME: Influenza A (H3N2) virus

Source: Allantoic fluid of 10 days old embryonated eggs, inoculated with influenza A virus, strain A/Panama/2007/99.

Purification: Ultracentrifugation with 10 – 40 % sucrose gradient.

Presentation: STE, 0.1 % sodium azide (NaN₃) and 0.005 % thimerosal.

Inactivation: Thimerosal and beta propiolactone treatment

This product has been treated in a manner consistent with methods of inactivation. Generally accepted good laboratory practices appropriate to microbiological/viral safe handling practices and techniques are required at work.

Analysis: Purity > 90 %

Immunological activity: Tested with anti-influenza A monoclonal antibodies in ELISA.

Application: Serological studies of influenza A virus, immunogen for antibody production.

Storage: - 20 °C

Material safety note: Contains sodium azide (0.1 %) as preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.