



**CATALOGUE #:** 8T24

**PRODUCT NAME:** Human skeletal muscle troponin T (skTnT)

---

**Source:** Human skeletal muscle.  
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

**Applications:** Immunogen for antibody production, troponin biochemical and immunochemical studies.

**Analysis:** Purity more than 95 % (SDS-PAGE).  
Protein content determined by the Lowry procedure.  
Immunological identity confirmed by reaction with monoclonal antibody that is specific for the TnT skeletal muscle isoform.

**Presentation:** Lyophilized from protein solution in 0.01 M HCl.

**Storage:** -20 °C (- 70 °C for long term storage)