



CATALOGUE #: 8T57

PRODUCT NAME: Human troponin C (TnC)

Source: Human heart tissue.
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

Applications: TnC is suitable for use as a standard, immunogen for antisera production, stabilizer of Tnl in solutions, troponin biochemical and immunochemical studies, tracer for iodination.

Analysis: Purity > 98 % (SDS-PAGE)

Presentation: Lyophilized from protein solution in 150 mM KCl, 20mM Tris-HCl, pH 7.5.
Product is recommended to reconstitute with distilled water.

Storage: -20 °C (- 70 °C for long term storage)