



CATALOGUE #: 8TG52

PRODUCT NAME: Human thyroglobulin

- Source:** Human thyroid gland.
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.
- Applications:** Human thyroglobulin is suitable for use as a standard in immunoassay for determination of serum thyroglobulin levels, immunogen for antisera production, and tracer for iodination.
- Analysis:** Purity > 90 % (SDS-PAGE).
Immunological identity confirmed by reaction with monoclonal antibody that is specific for the human thyroglobulin.
- Presentation:** Lyophilized from 25 mM Tris-HCl, pH8.
It is recommended to reconstitute the product in distilled water.
- Storage:** - 20 °C