

Data Sheet

D-Biotin

Cat. No.: 2-1016-002, 2-1016-005, 2-1016-010,
 2-1016-025

Lot No.: 1016 -

Description	Vitamin H; (+)-Biotin D-Biotin is an efficient competitor for the mild elution of <i>Strep</i> -tag II fusion proteins from <i>Strep</i> -Tactin columns. pH 8.0 is recommended. At least the buffer should be slightly alkaline for efficient dissolution. Important note: After elution of <i>Strep</i>-tag[®] proteins using D-Biotin the <i>Strep</i>-Tactin[®] resin can not be regenerated by <i>Strep</i>-tag Buffer R or HABA.
Form	white powder
FW	244.3, C ₁₀ H ₁₆ N ₂ O ₃ S
Stability	12 months after shipping
Storage	2 - 8 °C
Shipment	RT

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.